

Product Datasheet

Anti-Influenza A virus (H1N1)
V1301

Product Name	Anti-Influenza A virus (H1N1)
Catalog Number	V1301
Description	Purified IgG fraction prepared by methods developed by Medix Biochemica
Species/Host	Goat
Tested Applications	WB, ELISA, IFA, IHA, FACS
Brand	ViroStat
Form/Appearance	Purified IgG fraction prepared by methods developed by Medix Biochemica
Concentration	5 mg/ml
Storage	Short term (<6 months): +2-8 °C Long term: -20 °C
Isotype	IgG
Clonality	Polyclonal
Purity	≥ 95 %
Buffer	0.01M Phosphate Buffered Saline, pH 7.2, containing 0.1% sodium azide
Cross Reactivity	Non-reactive with influenza B, RSV, Para 1-3, and Adenovirus. Does not react with HEp-2 cells. May react with chicken cellular proteins
Specificity	Purified virions, H1N1 by IHA
Shelf Life	6 years