

Product Datasheet

Anti-Influenza A virus (H1N1) Biotin V1307

Product Name	Anti-Influenza A virus (H1N1) Biotin
Catalog Number	V1307
Species/Host	Goat
Tested Applications	ELISA, WB, IFA, IHA, FACS
Brand	ViroStat
Form/Appearance	Purified from antiserum and covalently coupled with the N-Hydroxysuccinimide ester of biotin under mild conditions to give a high degree of substitution.
Concentration	5 mg/ml
Storage	Short term (<6 months): 2-8°C Long term: -20°C
Isotype	IgG
Clonality	Polyclonal
Purity	≥95%
Buffer	0.01M Phosphate Buffered Saline, pH 7.2, containing 0.1% sodium azide
Cross Reactivity	Non-reactive with influenza B, RSV, Para 1-3, and Adenovirus. Does not react with Hep-2 cells. May react with chicken cellular proteins
Specificity	Purified virions, H1N1 by IHA
Shelf Life	Two years from date of manufacture