

Product Datasheet

Anti-Influenza A virus (H3N2) V1311

Product Name	Anti-Influenza A virus (H3N2)
Catalog Number	V1311
Species/Host	Goat
Tested Applications	Indirect immunofluorescence, IHA
Brand	ViroStat
Form/Appearance	Purified IgG fraction prepared by methods developed by Medix Biochemica
Concentration	5 mg/ml
Storage	Short term (<6 months): 2-8°C Long term: -20°C
Isotype	IgG
Clonality	Polyclonal
Purity	≥95%
Buffer	0.01M Phosphate Buffered Saline, pH 7.2, containing 0.1% sodium azide
Cross Reactivity	Non-reactive with Influenza B, RSV, Para 1-3, Adenovirus, & HEp-2 cells. May react with chicken cellular proteins
Specificity	Purified virions, H3N2 by IHA
Shelf Life	Six years from date of manufacture