

Product Datasheet

Anti-Influenza A virus (H3N2) FITC V1313

Product Name	Anti-Influenza A virus (H3N2) FITC
Catalog Number	V1313
Description	Purified from antiserum by protein A and covalently coupled with high purity Isomer I of fluorescein isothiocyanate (FITC)
Species/Host	Goat
Tested Applications	IFA, Indirect IFA, IHA
Brand	ViroStat
Form/Appearance	Purified from antiserum by protein A and covalently coupled with high purity Isomer I of fluorescein isothiocyanate (FITC)
Concentration	5 mg/ml
Storage	Short term (<6 months): +2-8 °C Long term: -20 °C
Isotype	IgG
Clonality	Polyclonal
Purity	≥ 95 %
Buffer	0.01M Phosphate Buffered Saline, pH 7.2, containing 0.1% sodium azide
Cross Reactivity	Non-reactive with Influenza B, RSV, Para 1-3, Adenovirus, & HEp-2 cells. May react with chicken cellular proteins
Specificity	Purified virions, H3N2 by IHA
Shelf Life	2 years