

## Product Datasheet

### Anti-Human IL-23 Monoclonal Antibody, Azide Free Clone B-F43 853.930.005

|                            |   |
|----------------------------|---|
| <b>Product Name</b>        | Anti-Human IL-23 Monoclonal Antibody, Azide Free Clone B-F43                                  |
| <b>Catalog Number</b>      | 853.930.005   |
| <b>Description</b>         | Mouse monoclonal  |
| <b>Species/Host</b>        | Mouse   |
| <b>Conjugation</b>         | Azide Free  |
| <b>Tested Applications</b> | ELISA   |
| <b>Brand</b>               | Diaclone  |
| <b>Form/Appearance</b>     | Phosphate-buffered saline. Sterile-filtered through 0.22 µm.<br>Carrier and preservative free |
| <b>Storage</b>             | Stable at +2-8°C for 12 months. For longer storage freeze aliquots.                           |
| <b>Isotype</b>             | Mouse IgG1 Kappa light chain  |
| <b>Clonality</b>           | Monoclonal  |
| <b>Source</b>              | Recombinant antibody from the original clone  |
| <b>Specificity</b>         | Recognizes both recombinant and natural human IL-23, p19 sub-unit                             |
| <b>Immunisation</b>        | Recombinant human IL-23 (sub-unit p19)  |