

# Product Datasheet

## Anti-Influenza A virus (matrix) V1321

Product Name	Anti-Influenza A virus (matrix)
Catalog Number	V1321
Species/Host	Goat
Tested Applications	Indirect immunoflourescence
Brand	ViroStat
Form/Appearance	Purified IgG fraction prepared by methods developed by Medix Biochemica
Concentration	5 mg/ml
Storage	Short term (<6 months): 2-8°C Long term: -20°C
Isotype	IgG
Clonality	Polyclonal
Purity	≥95%
Buffer	0.01M Phosphate Buffered Saline, pH 7.2, containing 0.1% sodium azide
Cross Reactivity	Non-reactive with the M2 matrix protein. Non-reactive with Hep-2 cells by indirect immunofluorescence
Specificity	Purified matrix protein (M1)
Shelf Life	Six years from date of manufacture