

## Product Datasheet

### Anti-Adenovirus FITC V1403

<b>Product Name</b>	Anti-Adenovirus FITC
<b>Catalog Number</b>	V1403
<b>Species/Host</b>	Goat
<b>Tested Applications</b>	ELISA, IFA, WB immunoassay
<b>Brand</b>	ViroStat
<b>Form/Appearance</b>	Purified from antiserum by protein A and covalently coupled with high purity Isomer I of fluorescein isothiocyanate (FITC)
<b>Concentration</b>	5 mg/ml
<b>Storage</b>	Short term (<6 months): 2-8°C Long term: -20°C
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Purity</b>	≥95%
<b>Buffer</b>	0.01M Phosphate Buffered Saline, pH 7.2, containing 0.1% sodium azide
<b>Cross Reactivity</b>	Reactive with Adeno types 1, 2, 3, 5, 6, 7a, 8, 31, 40, & 41. Others not tested. Non-reactive with Para 1-3, influenza A & B, RSV, hMPV, & Rhinovirus (type 16). Non-reactive with Hep-2, LLCMK2, & Vero cells by indirect immunofluorescence
<b>Specificity</b>	Numerous Adeno proteins including hexon
<b>Shelf Life</b>	Two years from date of manufacture