

Product Datasheet

Anti-Adenovirus HRP
V1404

Product Name	Anti-Adenovirus HRP
Catalog Number	V1404
Description	Purified from antiserum and covalently coupled to a highly purified preparation of horseradish peroxidase (HRP; RZ>3) by methods developed by Medix Biochemica
Species/Host	Goat
Tested Applications	ELISA
Brand	ViroStat
Form/Appearance	Purified from antiserum and covalently coupled to a highly purified preparation of horseradish peroxidase (HRP; RZ>3) by methods developed by Medix Biochemica
Concentration	1 mg/ml
Storage	Short term (<6 months): +2-8 °C Long term: -20 °C
Isotype	IgG
Clonality	Polyclonal
Purity	≥ 95 %
Buffer	0.01M Phosphate Buffered Saline, pH 7.2, containing bovine serum albumin (10 mg/ml) as a stabilizer and preserved with 0.002% thimerosal
Cross Reactivity	Reactive with Adeno types 1, 2, 3, 5, 6, 7a, 8, 31, 40, & 41. Others not tested. Non-reactive with Para 1-3, influenza A & B, RSV, hMPV, Rhinovirus (type 16), Hep-2, LLCMK2, & Vero cells
Specificity	Numerous Adeno proteins including hexon
Shelf Life	2 years