

Product Datasheet

Anti-Streptococcus Group A V1501

Product Name	Anti-Streptococcus Group A
Catalog Number	V1501
Description	Immunoaffinity purified specific antibody
Species/Host	Goat
Tested Applications	ELISA, LF, IFA
Brand	ViroStat
Form/Appearance	Immunoaffinity purified specific antibody
Concentration	≥ 2.0 mg/ml
Storage	Short term (<6 months): +2-8 °C Long term: -20 °C
Isotype	IgG
Clonality	Polyclonal
Purity	≥ 95 %
Buffer	0.01M Phosphate Buffered Saline, pH 7.2, containing 0.1% sodium azide
Cross Reactivity	Non-reactive with other Strep groups
Specificity	Type specific carbohydrate of group A Streptococcus
Shelf Life	6 years