

Product Datasheet

Anti-Human IL-17A Monoclonal Antibody, FITC Conjugated Clone B-A50 853.941.019

Product Name	Anti-Human IL-17A Monoclonal Antibody, FITC Conjugated Clone B-A50
Catalog Number	853.941.019
Description	Mouse monoclonal
Species/Host	Mouse
Conjugation	FITC
Tested Applications	FC
Immunogen	Recombinant human IL-17A
Brand	Diaclone
Form/Appearance	Phosphate-buffered saline with 5 % BSA and 0.1 % sodium azide
Storage	Stable at +2-8°C for 12 months. For longer storage freeze aliquots.
Note	Use 10 µl to label 10 ⁶ cells or 100 µl of whole blood
Isotype	Mouse IgG2b Kappa light chain
Clonality	Monoclonal
Specificity	Recognizes both natural and recombinant IL-17A
Hybridoma	Myeloma X63/AG.8653 x Balb/c node cells