

Product Datasheet

Anti-Streptococcus Group B FITC V1523

Product Name	Anti-Streptococcus Group B FITC
Catalog Number	V1523
Description	Purified from antiserum by protein A and covalently coupled with high purity Isomer I of fluorescein isothiocyanate (FITC)
Species/Host	Rabbit
Tested Applications	ELISA, LF, IFA
Brand	ViroStat
Form/Appearance	Purified from antiserum by protein A and covalently coupled with high purity Isomer I of fluorescein isothiocyanate (FITC)
Concentration	1 mg/ml
Storage	Short term (<6 months): +2-8 °C Long term: -20 °C
Isotype	IgG
Clonality	Polyclonal
Purity	≥ 95 %
Buffer	0.01M Phosphate Buffered Saline, pH 7.2, containing 0.1% sodium azide
Cross Reactivity	Non-reactive with other Strep groups (A, C, D, F, & G)
Specificity	Type specific carbohydrate of group B Streptococcus
Shelf Life	2 years