

# Product Datasheet

## Anti-Streptococcus Group B Biotin V1527

Product Name	Anti-Streptococcus Group B Biotin
Catalog Number	V1527
Species/Host	Rabbit
Tested Applications	LF, ELISA, IFA
Brand	ViroStat
Form/Appearance	Purified from antiserum and covalently coupled with the N-Hydroxysuccinimide ester of biotin under mild conditions to give a high degree of substitution.
Concentration	1 mg/ml
Storage	Short term (<6 months): 2-8°C Long term: -20°C
Isotype	IgG
Clonality	Polyclonal
Purity	≥95%
Buffer	0.01M Phosphate Buffered Saline, pH 7.2, containing 0.1% sodium azide
Cross Reactivity	Non-reactive with other Strep groups (A, C, D, F, & G)
Specificity	Type specific carbohydrate of group B Streptococcus
Shelf Life	Two years from date of manufacture