

## Product Datasheet

### Anti-Human IL-17A Monoclonal Antibody, PE Conjugated Clone B-A50 853.942.019

<b>Product Name</b>	Anti-Human IL-17A Monoclonal Antibody, PE Conjugated Clone B-A50
<b>Catalog Number</b>	853.942.019
<b>Description</b>	Mouse monoclonal
<b>Species/Host</b>	Mouse
<b>Conjugation</b>	PE
<b>Tested Applications</b>	FC
<b>Immunogen</b>	Recombinant human IL-17A
<b>Brand</b>	Diaclone
<b>Form/Appearance</b>	Phosphate-buffered saline with 5% BSA and 0.1% sodium azide
<b>Storage</b>	Stable at +2-8°C for six months after reconstitution. DO NOT FREEZE.
<b>Note</b>	Reconstitute with 1 ml deionised water Use 10 µl to label 10 <sup>6</sup> cells or 100 µl of whole blood
<b>Isotype</b>	Mouse IgG2b Kappa light chain
<b>Clonality</b>	Monoclonal
<b>Specificity</b>	Recognizes both natural and recombinant IL-17A
<b>Hybridoma</b>	Myeloma X63/AG.8653 x Balb/c node cells