

# Product Datasheet

## Anti-Chlamydia trachomatis V1601

Product Name	Anti-Chlamydia trachomatis
Catalog Number	V1601
Description	Purified IgG fraction prepared by methods developed by Medix Biochemica
Species/Host	Goat
Tested Applications	Indirect IFA
Brand	ViroStat
Form/Appearance	Purified IgG fraction prepared by methods developed by Medix Biochemica
Concentration	5 mg/ml
Storage	Short term (<6 months): +2-8 °C Long term: -20 °C
Isotype	IgG
Clonality	Polyclonal
Purity	≥ 95 %
Buffer	0.01M Phosphate Buffered Saline, pH 7.2, containing 0.1% sodium azide
Cross Reactivity	Reactive with Chlamydia psittacii & Chlamydia pneumoniae (TWAR). Non-reactive with HEp-2 cells and egg yolk sac
Specificity	Purified elementary bodies, disrupted
Shelf Life	6 years