

Product Datasheet

Anti-Chlamydia trachomatis FITC V1613

Product Name	Anti-Chlamydia trachomatis FITC
Catalog Number	V1613
Description	Purified from antiserum by protein A and covalently coupled with high purity Isomer I of fluorescein isothiocyanate (FITC)
Species/Host	Rabbit
Tested Applications	IFA, Indirect IFA
Brand	ViroStat
Form/Appearance	Purified from antiserum by protein A and covalently coupled with high purity Isomer I of fluorescein isothiocyanate (FITC)
Concentration	5 mg/ml
Storage	Short term (<6 months): +2-8 °C Long term: -20 °C
Isotype	IgG
Clonality	Polyclonal
Purity	≥ 95 %
Buffer	0.01M Phosphate Buffered Saline, pH 7.2, containing 0.1% sodium azide
Cross Reactivity	Reactive with Chlamydia psittacii & Chlamydia pneumoniae (TWAR). Non-reactive with Hep-2 cells and egg yolk sac
Specificity	Purified elementary bodies, disrupted
Shelf Life	2 years