

Product Datasheet

Anti-Chlamydia trachomatis (MOMP) HRP V1624

Product Name	Anti-Chlamydia trachomatis (MOMP) HRP
Catalog Number	V1624
Species/Host	Goat
Tested Applications	ELISA
Brand	ViroStat
Form/Appearance	Purified from antiserum and covalently coupled to a highly purified preparation of horseradish peroxidase (HRP; RZ>3) by methods developed by Medix Biochemica
Concentration	1 mg/ml
Storage	Short term (<6 months): 2-8°C Long term: -20°C
Isotype	IgG
Clonality	Polyclonal
Purity	≥95%
Buffer	0.01M Phosphate Buffered Saline, pH 7.2, containing bovine serum albumin (10 mg/ml) as a stabilizer and preserved with 0.002% thimerosal
Cross Reactivity	Non-reactive with C psittacii and C. pneumoniae in MIF. Negative vs HEp-2 cells and egg yolk sac
Specificity	Major Outer Membrane Protein (MOMP)
Shelf Life	Two years from date of manufacture