

Product Datasheet

Anti-Chlamydia trachomatis (MOMP) Biotin V1627

Product Name	Anti-Chlamydia trachomatis (MOMP) Biotin
Catalog Number	V1627
Species/Host	Goat
Tested Applications	ELISA, Indirect immunofluorescence
Brand	ViroStat
Form/Appearance	Purified from antiserum and covalently coupled with the N-Hydroxysuccinimide ester of biotin under mild conditions to give a high degree of substitution.
Concentration	5 mg/ml
Storage	Short term (<6 months): 2-8°C Long term: -20°C
Isotype	IgG
Clonality	Polyclonal
Purity	≥95%
Buffer	0.01M Phosphate Buffered Saline, pH 7.2, containing 0.1% sodium azide
Cross Reactivity	Non-reactive with C psittacii and C. pneumoniae in MIF. Negative vs HEp-2 cells and egg yolk sac
Specificity	Major Outer Membrane Protein (MOMP)
Shelf Life	Two years from date of manufacture