

# Product Datasheet

## Anti-HIV-1 V1901

|                     |   |
|---------------------|---|
| Product Name        | Anti-HIV-1  |
| Catalog Number      | V1901   |
| Species/Host        | Goat  |
| Tested Applications | Indirect immunofluorescence   |
| Brand               | ViroStat  |
| Form/Appearance     | Purified IgG fraction prepared by methods developed by Medix Biochemica |
| Concentration       | 5 mg/ml   |
| Storage             | Short term (<6 months): 2-8°C Long term: -20°C                          |
| Isotype             | IgG   |
| Clonality           | Polyclonal  |
| Purity              | ≥95%  |
| Buffer              | 0.01M Phosphate Buffered Saline, pH 7.2, containing 0.1% sodium azide   |
| Cross Reactivity    | Non-reactive with T4+ human lymphocytes                                 |
| Specificity         | Purified virions, disrupted   |
| Shelf Life          | Six years from date of manufacture                                      |