

Product Datasheet

Anti-HIV-1 (p24) FITC V1953

Product Name	Anti-HIV-1 (p24) FITC
Catalog Number	V1953
Species/Host	Goat
Tested Applications	IFA, WB immunoassay
Brand	ViroStat
Form/Appearance	Purified from antiserum by protein A and covalently coupled with high purity Isomer I of fluorescein isothiocyanate (FITC)
Concentration	5 mg/ml
Storage	Short term (<6 months): 2-8°C Long term: -20°C
Isotype	IgG
Clonality	Polyclonal
Purity	≥95%
Buffer	0.01M Phosphate Buffered Saline, pH 7.2, containing 0.1% sodium azide
Cross Reactivity	Non-reactive with human T or B cells. Non-reactive with human or bovine sera
Specificity	p24 (predominate p24 band pattern in WB)
Shelf Life	Two years from date of manufacture