

Product Datasheet

Anti-HIV-1 (p24) Biotin V1957

Product Name	Anti-HIV-1 (p24) Biotin
Catalog Number	V1957
Description	Purified from antiserum and covalently coupled with the N-Hydroxysuccinimide ester of biotin under mild conditions to give a high degree of substitution.
Species/Host	Goat
Tested Applications	WB, ELISA
Brand	ViroStat
Form/Appearance	Purified from antiserum and covalently coupled with the N-Hydroxysuccinimide ester of biotin under mild conditions to give a high degree of substitution.
Concentration	5 mg/ml
Storage	Short term (<6 months): +2-8 °C Long term: -20 °C
Isotype	IgG
Clonality	Polyclonal
Purity	≥ 95 %
Buffer	0.01M Phosphate Buffered Saline, pH 7.2, containing 0.1% sodium azide
Cross Reactivity	Non-reactive with human T or B cells. Non-reactive with human or bovine sera
Specificity	p24 (predominate p24 band pattern in WB)
Shelf Life	2 years