

Product Datasheet

Anti-HIV-1 (gp41) V1971

| | |
|---------------------|--|
| Product Name | Anti-HIV-1 (gp41) |
| Catalog Number | V1971 |
| Description | Purified IgG fraction prepared by methods developed by Medix Biochemica |
| Species/Host | Goat |
| Tested Applications | WB |
| Brand | ViroStat |
| Form/Appearance | Purified IgG fraction prepared by methods developed by Medix Biochemica |
| Concentration | 5 mg/ml |
| Storage | Short term (<6 months): +2-8 °C Long term: -20 °C |
| Isotype | IgG |
| Clonality | Polyclonal |
| Purity | ≥ 95 % |
| Buffer | 0.01M Phosphate Buffered Saline, pH 7.2, containing 0.1% sodium azide |
| Cross Reactivity | Non-reactive with human T or B cells. Reactive with gp41 in WB immunoassay |
| Specificity | a.a.546-a.a.682 sequence from strain HxB2 |
| Shelf Life | 6 years |