

Product Datasheet

Anti-Human CD213a2 Monoclonal Antibody, Azide Free Clone B-A34 853.960

Product Name	Anti-Human CD213a2 Monoclonal Antibody, Azide Free Clone B-A34
Catalog Number	853.960
Description	Mouse monoclonal antibody. Source: Myeloma X63/AG.8653 x Balb/c spleen cells
Species/Host	Mouse
Conjugation	Azide Free
Tested Applications	ELISA
Immunogen	IL-13Ra2 transfected CHO cell line
Alternative Names	IL-13Ra2
Brand	Diaclone
Form/Appearance	Liquid
Concentration	1.0 mg/ml (+/- 10 %)
Storage	Stable at +2-8 °C for 12 months. For longer storage freeze aliquots.
Isotype	Mouse IgG2a Kappa light chain
Clonality	Monoclonal
Buffer	Phosphate-buffered saline. Sterile-filtered through 0.22µm. Carrier and preservative free.
Specificity	Recognises the IL-13 receptor a2 (IL-13 Ra2), a 74 kDa protein
Hybridoma	Myeloma X63/AG.8653 x Balb/c spleen cells