

Product Datasheet

Anti-Bacillus FITC V3403

Product Name Anti-Bacillus FITC

Catalog Number V3403

Species/Host Rabbit

Tested Applications IFA

Brand ViroStat

Form/Appearance Purified from antiserum by protein A and covalently coupled with high purity Isomer I of fluorescein isothiocyanate (FITC)

Concentration 5 mg/ml

Storage Short term (<6 months): 2-8°C Long term: -20°C

Isotype IgG

Clonality Polyclonal

Purity ≥95%

Buffer 0.01M Phosphate Buffered Saline, pH 7.2, containing 0.1% sodium azide

Cross Reactivity Antiserum is unabsorbed and is expected to cross react with various Bacillus species

Specificity Various Bacillus species

Shelf Life Two years from date of manufacture