

# Product Datasheet

## Anti-Diphtheria toxin FITC V4703

Product Name	Anti-Diphtheria toxin FITC
Catalog Number	V4703
Species/Host	Goat
Tested Applications	IFA
Brand	ViroStat
Form/Appearance	Purified from antiserum by protein A and covalently coupled with high purity Isomer I of fluorescein isothiocyanate (FITC)
Concentration	5 mg/ml
Storage	Short term (<6 months): 2-8°C Long term: -20°C
Isotype	IgG
Clonality	Polyclonal
Purity	≥95%
Buffer	0.01M Phosphate Buffered Saline, pH 7.2, containing 0.1% sodium azide
Cross Reactivity	Demonstrated reactivity with purified toxin
Specificity	Purified diphtheria toxin/toxoid
Shelf Life	Two years from date of manufacture