

Product Datasheet

Anti-Legionella pneumophila FITC
V6053

Product Name	Anti-Legionella pneumophila FITC
Catalog Number	V6053
Species/Host	Rabbit
Tested Applications	LF, ELISA, IFA
Brand	ViroStat
Form/Appearance	Purified from antiserum by protein A and covalently coupled with high purity Isomer I of fluorescein isothiocyanate (FITC)
Concentration	5 mg/ml
Storage	Short term (<6 months): 2-8°C Long term: -20°C
Isotype	IgG
Clonality	Polyclonal
Purity	≥95%
Buffer	0.01M Phosphate Buffered Saline, pH 7.2, containing 0.1% sodium azide
Cross Reactivity	Antiserum is not absorbed and may cross react with related microorganisms. Reactive with L. pneumophila serogroups 1-12 in IFA
Specificity	All antigens of intact microorganism
Shelf Life	Two years from date of manufacture