

Product Datasheet

Anti-Legionella pneumophila Biotin V6057

Product Name	Anti-Legionella pneumophila Biotin
Catalog Number	V6057
Description	Purified from antiserum and covalently coupled with the N-Hydroxysuccinimide ester of biotin under mild conditions to give a high degree of substitution.
Species/Host	Rabbit
Tested Applications	ELISA, LF, IFA
Brand	ViroStat
Form/Appearance	Purified from antiserum and covalently coupled with the N-Hydroxysuccinimide ester of biotin under mild conditions to give a high degree of substitution.
Concentration	5 mg/ml
Storage	Short term (<6 months): +2-8 °C Long term: -20 °C
Isotype	IgG
Clonality	Polyclonal
Purity	≥ 95 %
Buffer	0.01M Phosphate Buffered Saline, pH 7.2, containing 0.1% sodium azide
Cross Reactivity	Antiserum is not absorbed and may cross react with related microorganisms. Reactive with L. pneumophila serogroups 1-12 in IFA
Specificity	All antigens of intact microorganism
Shelf Life	2 years