

## Product Datasheet

### Anti-Human IL-17F Monoclonal Antibody, Azide Free Clone B-G46 853.990.005

<b>Product Name</b>	Anti-Human IL-17F Monoclonal Antibody, Azide Free Clone B-G46
<b>Catalog Number</b>	853.990.005
<b>Description</b>	Mouse monoclonal
<b>Species/Host</b>	Mouse
<b>Conjugation</b>	Azide Free
<b>Tested Applications</b>	ELISA, ELISpot
<b>Brand</b>	Diaclone
<b>Form/Appearance</b>	Phosphate-buffered saline. Sterile-filtered through 0.22 µm. Carrier and preservative free
<b>Storage</b>	Stable at +2-8°C for 12 months. For longer storage freeze aliquots.
<b>Isotype</b>	Mouse IgG2b Kappa light chain
<b>Clonality</b>	Monoclonal
<b>Source</b>	Recombinant human IL-17F
<b>Specificity</b>	Recognizes both recombinant and natural human IL-17F
<b>Hybridoma</b>	Myeloma X63/AG.8653 x Balb/c node cells
<b>Immunisation</b>	Recombinant human IL-17F