

Product Datasheet

Anti-Neisseria meningitidis V6121

Product Name	Anti-Neisseria meningitidis
Catalog Number	V6121
Description	Purified from antiserum by protein A chromatography
Species/Host	Rabbit
Tested Applications	Indirect IFA
Brand	ViroStat
Form/Appearance	Purified from antiserum by protein A chromatography
Concentration	5 mg/ml
Storage	Short term (<6 months): +2-8 °C Long term: -20 °C
Isotype	IgG
Clonality	Polyclonal
Purity	≥ 95 %
Buffer	0.01M Phosphate Buffered Saline, pH 7.2, containing 0.1% sodium azide
Cross Reactivity	Reactive with serogroups A, B, & C. Other serogroups not tested. Antiserum has not been absorbed and may react with related microorganisms
Specificity	All antigens
Shelf Life	6 years