

Product Datasheet

Anti-Neisseria meningitidis Biotin V6127

Product Name	Anti-Neisseria meningitidis Biotin
Catalog Number	V6127
Species/Host	Rabbit
Tested Applications	ELISA, Indirect IFA
Brand	ViroStat
Form/Appearance	Purified from antiserum and covalently coupled with the N-Hydroxysuccinimide ester of biotin under mild conditions to give a high degree of substitution.
Concentration	5 mg/ml
Storage	Short term (<6 months): 2-8°C Long term: -20°C
Isotype	IgG
Clonality	Polyclonal
Purity	≥95%
Buffer	0.01M Phosphate Buffered Saline, pH 7.2, containing 0.1% sodium azide
Cross Reactivity	Reactive with serogroups A, B, & C. Other serogroups not tested. Antiserum has not been absorbed and may react with related microorganisms
Specificity	All antigens
Shelf Life	Two years from date of manufacture