

## Product Datasheet

### Anti-Candida albicans V6411

<b>Product Name</b>	Anti-Candida albicans
<b>Catalog Number</b>	V6411
<b>Description</b>	Purified from antiserum by protein A chromatography
<b>Species/Host</b>	Rabbit
<b>Tested Applications</b>	IHC, ELISA, IFA, Double-diffusion, CIE
<b>Brand</b>	ViroStat
<b>Form/Appearance</b>	Purified from antiserum by protein A chromatography
<b>Concentration</b>	5 mg/ml
<b>Storage</b>	Short term (<6 months): +2-8 °C Long term: -20 °C
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Purity</b>	≥ 95 %
<b>Buffer</b>	0.01M Phosphate Buffered Saline, pH 7.2, containing 0.1% sodium azide
<b>Cross Reactivity</b>	Non-reactive with human serum, urine, & spinal fluid. Antiserum has not been absorbed and does cross react with other yeasts
<b>Specificity</b>	Recognizes numerous proteins in a soluble C. albicans extract (IEP)
<b>Shelf Life</b>	6 years