

Product Datasheet

Anti-Candida albicans HRP V6414

Product Name	Anti-Candida albicans HRP
Catalog Number	V6414
Description	Purified from antiserum and covalently coupled to a highly purified preparation of horseradish peroxidase (HRP; RZ>3) by methods developed by Medix Biochemica
Species/Host	Rabbit
Tested Applications	ELISA
Brand	ViroStat
Form/Appearance	Purified from antiserum and covalently coupled to a highly purified preparation of horseradish peroxidase (HRP; RZ>3) by methods developed by Medix Biochemica
Concentration	1 mg/ml
Storage	Short term (<6 months): +2-8 °C Long term: -20 °C
Isotype	IgG
Clonality	Polyclonal
Purity	≥ 95 %
Buffer	0.01M Phosphate Buffered Saline, pH 7.2, containing bovine serum albumin (10 mg/ml) as a stabilizer and preserved with 0.002% thimerosal
Cross Reactivity	Non-reactive with human serum, urine, & spinal fluid. Antiserum has not been absorbed and does cross react with other yeasts
Specificity	Recognizes numerous proteins in a soluble C. albicans extract (IEP)
Shelf Life	2 years