

Product Datasheet

Anti-Candida albicans Biotin V6417

Product Name	Anti-Candida albicans Biotin
Catalog Number	V6417
Description	Purified from antiserum and covalently coupled with the N-Hydroxysuccinimide ester of biotin under mild conditions to give a high degree of substitution.
Species/Host	Rabbit
Tested Applications	IHC, ELISA, IFA, Double-diffusion, CIE
Brand	ViroStat
Form/Appearance	Purified from antiserum and covalently coupled with the N-Hydroxysuccinimide ester of biotin under mild conditions to give a high degree of substitution.
Concentration	5 mg/ml
Storage	Short term (<6 months): +2-8 °C Long term: -20 °C
Isotype	IgG
Clonality	Polyclonal
Purity	≥ 95 %
Buffer	0.01M Phosphate Buffered Saline, pH 7.2, containing 0.1% sodium azide
Cross Reactivity	Non-reactive with human serum, urine, & spinal fluid. Antiserum has not been absorbed and does cross react with other yeasts
Specificity	Recognizes numerous proteins in a soluble C. albicans extract (IEP)
Shelf Life	2 years