

# Product Datasheet

## Anti-Helicobacter pylori HRP V6604

|                     |  |
|---------------------|--|
| Product Name        | Anti-Helicobacter pylori HRP   |
| Catalog Number      | V6604  |
| Species/Host        | Rabbit   |
| Tested Applications | ELISA  |
| Brand               | ViroStat   |
| Form/Appearance     | Purified from antiserum and covalently coupled to a highly purified preparation of horseradish peroxidase (HRP; RZ>3) by methods developed by Medix Biochemica |
| Concentration       | 1 mg/ml  |
| Storage             | Short term (<6 months): 2-8°C Long term: -20°C   |
| Isotype             | IgG  |
| Clonality           | Polyclonal   |
| Purity              | ≥95%   |
| Buffer              | 0.01M Phosphate Buffered Saline, pH 7.2, containing bovine serum albumin (10 mg/ml) as a stabilizer and preserved with 0.002% thimerosal                       |
| Cross Reactivity    | Antiserum has not been absorbed and may react with related microorganisms  |
| Specificity         | Whole cell lysate  |
| Shelf Life          | Two years from date of manufacture   |