

Product Datasheet

Anti-Aspergillus V6821

Product Name	Anti-Aspergillus
Catalog Number	V6821
Description	Purified from antiserum by protein A chromatography
Species/Host	Rabbit
Tested Applications	IHC
Brand	ViroStat
Form/Appearance	Purified from antiserum by protein A chromatography
Concentration	5 mg/ml
Storage	Short term (<6 months): +2-8 °C Long term: -20 °C
Isotype	IgG
Clonality	Polyclonal
Purity	≥ 95 %
Buffer	0.01M Phosphate Buffered Saline, pH 7.2, containing 0.1% sodium azide
Cross Reactivity	Reactive with soluble proteins from common Aspergillus species. Antiserum is not absorbed and may react with related microorganisms
Specificity	Soluble proteins from common Aspergillus species
Shelf Life	6 years