

Product Datasheet

Anti-Aspergillus FITC
V6823

Product Name	Anti-Aspergillus FITC
Catalog Number	V6823
Description	Purified from antiserum by protein A and covalently coupled with high purity Isomer I of fluorescein isothiocyanate (FITC)
Species/Host	Rabbit
Tested Applications	IHC, IFA
Brand	ViroStat
Form/Appearance	Purified from antiserum by protein A and covalently coupled with high purity Isomer I of fluorescein isothiocyanate (FITC)
Concentration	5 mg/ml
Storage	Short term (<6 months): +2-8 °C Long term: -20 °C
Isotype	IgG
Clonality	Polyclonal
Purity	≥ 95 %
Buffer	0.01M Phosphate Buffered Saline, pH 7.2, containing 0.1% sodium azide
Cross Reactivity	Reactive with soluble proteins from common Aspergillus species. Antiserum is not absorbed and may react with related microorganisms
Specificity	Soluble proteins from common Aspergillus species
Shelf Life	2 years