

# Product Datasheet

## Anti-Aspergillus HRP V6824

Product Name	Anti-Aspergillus HRP
Catalog Number	V6824
Species/Host	Rabbit
Tested Applications	ELISA
Brand	ViroStat
Form/Appearance	Purified from antiserum and covalently coupled to a highly purified preparation of horseradish peroxidase (HRP; RZ>3) by methods developed by Medix Biochemica
Concentration	1 mg/ml
Storage	Short term (<6 months): 2-8°C Long term: -20°C
Isotype	IgG
Clonality	Polyclonal
Purity	≥95%
Buffer	0.01M Phosphate Buffered Saline, pH 7.2, containing bovine serum albumin (10 mg/ml) as a stabilizer and preserved with 0.002% thimerosal
Cross Reactivity	Reactive with soluble proteins from common Aspergillus species. Antiserum is not absorbed and may react with related microorganisms
Specificity	Soluble proteins from common Aspergillus species
Shelf Life	Two years from date of manufacture