

Product Datasheet

Anti-Aspergillus Biotin V6827

Product Name	Anti-Aspergillus Biotin
Catalog Number	V6827
Description	Purified from antiserum and covalently coupled with the N-Hydroxysuccinimide ester of biotin under mild conditions to give a high degree of substitution.
Species/Host	Rabbit
Tested Applications	IHC, ELISA
Brand	ViroStat
Form/Appearance	Purified from antiserum and covalently coupled with the N-Hydroxysuccinimide ester of biotin under mild conditions to give a high degree of substitution.
Concentration	5 mg/ml
Storage	Short term (<6 months): +2-8 °C Long term: -20 °C
Isotype	IgG
Clonality	Polyclonal
Purity	≥ 95 %
Buffer	0.01M Phosphate Buffered Saline, pH 7.2, containing 0.1% sodium azide
Cross Reactivity	Reactive with soluble proteins from common Aspergillus species. Antiserum is not absorbed and may react with related microorganisms
Specificity	Soluble proteins from common Aspergillus species
Shelf Life	2 years