

Product Datasheet

Anti-Clostridia V6871

| | |
|---------------------|---|
| Product Name | Anti-Clostridia |
| Catalog Number | V6871 |
| Species/Host | Rabbit |
| Tested Applications | IFA |
| Brand | ViroStat |
| Form/Appearance | Purified from antiserum and covalently coupled with the N-Hydroxysuccinimide ester of biotin under mild conditions to give a high degree of substitution. |
| Concentration | 5 mg/ml |
| Storage | Short term (<6 months): 2-8°C Long term: -20°C |
| Isotype | IgG |
| Clonality | Polyclonal |
| Purity | ≥95% |
| Buffer | 0.01M Phosphate Buffered Saline, pH 7.2, containing 0.1% sodium azide |
| Cross Reactivity | Reactive with proteins from common Clostridia species. Reactive in IFA with Cl. Perfringens, Cl. Sordellii, Cl. Novyi, Cl. Septicum, Cl. Chauvoei |
| Specificity | Proteins from common Clostridia species |
| Shelf Life | Six years from date of manufacture |