

Product Datasheet

Anti-Clostridia FITC
V6873

Product Name	Anti-Clostridia FITC
Catalog Number	V6873
Description	Purified from antiserum by protein A and covalently coupled with high purity Isomer I of fluorescein isothiocyanate (FITC)
Species/Host	Rabbit
Tested Applications	IFA
Brand	ViroStat
Form/Appearance	Purified from antiserum by protein A and covalently coupled with high purity Isomer I of fluorescein isothiocyanate (FITC)
Concentration	5 mg/ml
Storage	Short term (<6 months): +2-8 °C Long term: -20 °C
Isotype	IgG
Clonality	Polyclonal
Purity	≥ 95 %
Buffer	0.01M Phosphate Buffered Saline, pH 7.2, containing 0.1% sodium azide
Cross Reactivity	Reactive with proteins from common Clostridia species. Reactive in IFA with Cl. Perfringens, Cl. Sordellii, Cl. Novyi, Cl. Septicum, Cl. Chauvoei
Specificity	Proteins from common Clostridia species
Shelf Life	2 years