

# Product Datasheet

## Anti-Clostridia HRP V6874

Product Name	Anti-Clostridia HRP
Catalog Number	V6874
Species/Host	Rabbit
Tested Applications	ELISA
Brand	ViroStat
Form/Appearance	Purified from antiserum and covalently coupled to a highly purified preparation of horseradish peroxidase (HRP; RZ>3) by methods developed by Medix Biochemica
Concentration	1 mg/ml
Storage	Short term (<6 months): 2-8°C Long term: -20°C
Isotype	IgG
Clonality	Polyclonal
Purity	≥95%
Buffer	0.01M Phosphate Buffered Saline, pH 7.2, containing bovine serum albumin (10 mg/ml) as a stabilizer and preserved with 0.002% thimerosal
Cross Reactivity	Reactive with proteins from common Clostridia species. Reactive with Cl. Perfringens, Cl. Sordellii, Cl. Novyi, Cl. Septicum, Cl. Chauvoei
Specificity	Proteins from common Clostridia species
Shelf Life	Two years from date of manufacture