

Product Datasheet

Anti-Clostridia Biotin V6877

Product Name	Anti-Clostridia Biotin
Catalog Number	V6877
Description	Purified from antiserum and covalently coupled with the N-Hydroxysuccinimide ester of biotin under mild conditions to give a high degree of substitution.
Species/Host	Rabbit
Tested Applications	ELISA, IFA
Brand	ViroStat
Form/Appearance	Purified from antiserum and covalently coupled with the N-Hydroxysuccinimide ester of biotin under mild conditions to give a high degree of substitution.
Concentration	5 mg/ml
Storage	Short term (<6 months): +2-8 °C Long term: -20 °C
Isotype	IgG
Clonality	Polyclonal
Purity	≥ 95 %
Buffer	0.01M Phosphate Buffered Saline, pH 7.2, containing 0.1% sodium azide
Cross Reactivity	Reactive with proteins from common Clostridia species. Reactive in IFA with Cl. Perfringens, Cl. Sordellii, Cl. Novyi, Cl. Septicum, Cl. Chauvoei
Specificity	Proteins from common Clostridia species
Shelf Life	2 years