

## Product Datasheet

### Anti-Staphylococcus aureus HRP V6884

<b>Product Name</b>	Anti-Staphylococcus aureus HRP
<b>Catalog Number</b>	V6884
<b>Description</b>	Purified from antiserum and covalently coupled to a highly purified preparation of horseradish peroxidase (HRP; RZ>3) by methods developed by Medix Biochemica
<b>Species/Host</b>	Rabbit
<b>Tested Applications</b>	ELISA
<b>Brand</b>	ViroStat
<b>Form/Appearance</b>	Purified from antiserum and covalently coupled to a highly purified preparation of horseradish peroxidase (HRP; RZ>3) by methods developed by Medix Biochemica
<b>Concentration</b>	1 mg/ml
<b>Storage</b>	Short term (<6 months): +2-8 °C Long term: -20 °C
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Purity</b>	≥ 95 %
<b>Buffer</b>	0.01M Phosphate Buffered Saline, pH 7.2, containing bovine serum albumin (10 mg/ml) as a stabilizer and preserved with 0.002% thimerosal
<b>Cross Reactivity</b>	Antiserum is not absorbed for related microorganisms
<b>Specificity</b>	Soluble and structural antigens of the whole bacterium
<b>Shelf Life</b>	2 years