

Product Datasheet

Anti-Klebsiella V6891

Product Name	Anti-Klebsiella
Catalog Number	V6891
Description	Purified from antiserum by protein A chromatography
Species/Host	Rabbit
Tested Applications	IFA
Brand	ViroStat
Form/Appearance	Purified from antiserum by protein A chromatography
Concentration	5 mg/ml
Storage	Short term (<6 months): +2-8 °C Long term: -20 °C
Isotype	IgG
Clonality	Polyclonal
Purity	≥ 95 %
Buffer	0.01M Phosphate Buffered Saline, pH 7.2, containing 0.1% sodium azide
Cross Reactivity	Antiserum is unabsorbed and may react with related Enterobacteriaceae
Specificity	Polyvalent for multiple Klebsiella antigens
Shelf Life	6 years