

## Product Datasheet

### Anti-Klebsiella FITC V6893

**Product Name** Anti-Klebsiella FITC

**Catalog Number** V6893

**Species/Host** Rabbit

**Tested Applications** IFA

**Brand** ViroStat

**Form/Appearance** Purified from antiserum by protein A and covalently coupled with high purity Isomer I of fluorescein isothiocyanate (FITC)

**Concentration** 5 mg/ml

**Storage** Short term (<6 months): 2-8°C Long term: -20°C

**Isotype** IgG

**Clonality** Polyclonal

**Purity** ≥95%

**Buffer** 0.01M Phosphate Buffered Saline, pH 7.2, containing 0.1% sodium azide

**Cross Reactivity** Antiserum is unabsorbed and may react with related Enterobacteriaceae

**Specificity** Polyvalent for multiple Klebsiella antigens

**Shelf Life** Two years from date of manufacture