

Product Datasheet

Anti-Lipid A V6921

Product Name	Anti-Lipid A
Catalog Number	V6921
Description	Purified IgG fraction prepared by methods developed by Medix Biochemica
Species/Host	Goat
Tested Applications	IFA
Brand	ViroStat
Form/Appearance	Purified IgG fraction prepared by methods developed by Medix Biochemica
Concentration	5 mg/ml
Storage	Short term (<6 months): +2-8 °C Long term: -20 °C
Isotype	IgG
Clonality	Polyclonal
Purity	≥ 95 %
Buffer	0.01M Phosphate Buffered Saline, pH 7.2, containing 0.1% sodium azide
Cross Reactivity	Reactive with Pseudomonas aeruginosa, Klebsiella pneumoniae, E. coli O157, Salmonella enteritidis, Enterobacter aerogenes, E. hermannii, Yersinia enterocolitica, Shigella sonnei
Specificity	LPS
Shelf Life	6 years