

Product Datasheet

Anti-Cryptosporidium parvum V7651

Product Name	Anti-Cryptosporidium parvum
Catalog Number	V7651
Description	Purified IgG fraction prepared by methods developed by Medix Biochemica
Species/Host	Goat
Tested Applications	ELISA, IFA, Indirect IFA
Brand	ViroStat
Form/Appears as	Purified IgG fraction prepared by methods developed by Medix Biochemica
Concentration	5 mg/ml
Storage	Short term (<6 months): +2-8 °C Long term: -20 °C
Isotype	IgG
Clonality	Polyclonal
Purity	≥ 95 %
Buffer	0.01M Phosphate Buffered Saline, pH 7.2, containing 0.1% sodium azide
Cross Reactivity	Reacts with a large rod-shaped microorganism found in bovine feces
Specificity	Intact Oocysts
Shelf Life	6 years