

Product Datasheet

Anti-Cholera toxin (Beta subunit) FITC V7923

Product Name	Anti-Cholera toxin (Beta subunit) FITC
Catalog Number	V7923
Description	Purified from antiserum by protein A and covalently coupled with high purity Isomer I of fluorescein isothiocyanate (FITC)
Species/Host	Rabbit
Tested Applications	ELISA, IFA
Brand	ViroStat
Form/Appearance	Purified from antiserum by protein A and covalently coupled with high purity Isomer I of fluorescein isothiocyanate (FITC)
Concentration	5 mg/ml
Storage	Short term (<6 months): +2-8 °C Long term: -20 °C
Isotype	IgG
Clonality	Polyclonal
Purity	≥ 95 %
Buffer	0.01M Phosphate Buffered Saline, pH 7.2, containing 0.1% sodium azide
Cross Reactivity	Reactive with E. coli heat labile toxin
Specificity	B subunit of Cholera toxin
Shelf Life	2 years