

Product Datasheet

Anti-Vaccinia virus V8101

Product Name	Anti-Vaccinia virus
Catalog Number	V8101
Description	Purified from antiserum by protein A chromatography
Species/Host	Rabbit
Tested Applications	IHC, IFA, PRNT, Indirect ELISA
Brand	ViroStat
Form/Appearance	Purified from antiserum by protein A chromatography
Concentration	5 mg/ml
Storage	Short term (<6 months): +2-8 °C Long term: -20 °C
Isotype	IgG
Clonality	Polyclonal
Purity	≥ 95 %
Buffer	0.01M Phosphate Buffered Saline, pH 7.2, containing 0.1% sodium azide
Cross Reactivity	Reactive with the Lister & MVA strains of Vaccinia and Monkeypox. Reactive with the B5R and A33R proteins of the virus. Non-reactive with Parainfluenza (1-3), RSV, Adeno, Influenza A&B, HS1 & uninfected cells
Specificity	Purified virions
Shelf Life	6 years