

Product Datasheet

Anti-Vaccinia virus HRP V8104

Product Name	Anti-Vaccinia virus HRP
Catalog Number	V8104
Description	Purified from antiserum and covalently coupled to a highly purified preparation of horseradish peroxidase (HRP; RZ>3) by methods developed by Medix Biochemica
Species/Host	Rabbit
Tested Applications	ELISA
Brand	ViroStat
Form/Appearance	Purified from antiserum and covalently coupled to a highly purified preparation of horseradish peroxidase (HRP; RZ>3) by methods developed by Medix Biochemica
Concentration	1 mg/ml
Storage	Short term (<6 months): +2-8 °C Long term: -20 °C
Isotype	IgG
Clonality	Polyclonal
Purity	≥ 95 %
Buffer	0.01M Phosphate Buffered Saline, pH 7.2, containing bovine serum albumin (10 mg/ml) as a stabilizer and preserved with 0.002% thimerosal
Cross Reactivity	Reactive with the Lister & MVA strains of Vaccinia and Monkeypox. Reactive with the B5R and A33R proteins of the virus. Non-reactive with Parainfluenza (1-3), RSV, Adeno, Influenza A&B, HS1 & uninfected cells
Specificity	Purified virions
Shelf Life	2 years