

Product Datasheet

Anti-Herpes simplex virus, type 2 (nuclear) V0115

Product Name	Anti-Herpes simplex virus, type 2 (nuclear)
Catalog Number	V0115
Tested Applications	IFA
Brand	ViroStat
Form/Appearance	Purified from mouse ascites by protein A chromatography
Concentration	≥ 1 mg/ml
Storage	2-8°C or -20°C. Avoid freeze-thaws. Maximum one freeze-thaw during storage.
Isotype	IgG2a
Clonality	Monoclonal
Purity	≥ 90%
Buffer	0.01M Phosphate Buffered Saline, pH 7.2, containing 0.1% sodium azide
Cross Reactivity	Reactive with HSV 2
Specificity	Nuclear (regulatory) protein
Shelf Life	Four years from date of manufacture